

ORAL PRESENTATION

Open Access

Prior epigenetic priming of cytokine genes in naive T cells is required for their subsequent activation by inducible enhancers

Sally R James¹, Sarah L Bevington², Fabio Mirabella¹, Marjorie Boissinot¹, Euan W Baxter¹, Sarion R Bowers¹, Peter N Cockerill^{1,2*}

From Epigenetics and Chromatin: Interactions and processes Boston, MA, USA. 11-13 March 2013

The inducible IL-3 and GM-CSF genes can only be efficiently expressed in T cells once they have been through a previous cycle of activation and undergone a process termed blast cell transformation [1]. An initial stimulus is required to bring naive T cells out of the resting state and into the cell cycle. This process is triggered by T cell receptor (TCR) stimulation, takes about 24 hours, and is associated with extensive nuclear remodeling. Once primed by a cycle of activation, T blast cells can maintain their ability to express IL-3 and GM-CSF for many cell divisions without the continual need for additional stimuli. In contrast, the IL-3 and GM-CSF genes cannot be induced in naive T cells that have never received a TCR stimulus. We show that this pattern of regulation of the IL-3/GM-CSF locus is controlled at two distinct levels:

- (1) During blast cell transformation, the IL-3/GM-CSF locus acquires an extensive array of DNase I Hypersensitive Sites (DHSs) which are then maintained indefinitely for many cell cycles [1,2]. These primed DHSs are marked by me2K4 histone H3, and they also persist in non-dividing memory T cells in the peripheral blood [1]. These DHSs are absent in the thymus, spleen T cells, and naive T cells in the blood. These DHSs do not function as classical enhancers, and we propose that they serve to maintain an active chromatin structure in previously activated T cells.
- (2) The expression of the lL-3 and GM-CSF genes is in each case dependent on the activation of inducible upstream enhancers. The lL-3 and GM-CSF enhancers appear as inducible DHSs T blast cells within 20 min of stimulation, and only acquire the me2K4H3 modification after stimulation. These enhancers are dependant on the

TCR inducible factors NFAT and AP-1. However, although AP-1 and NFAT family mRNAs are efficiently expressed in both naive T cells and in T blast cells, the enhancers only respond to induction of these factors in T blast cells. The inducible DHSs remain completely undetectable in naive T cells even after 4 hours of stimulation with direct activators of TCR signaling pathways (PMA and Calcium ionophore).

We propose that T blast cells and memory T cells have developed a strategy of maintaining a discrete class of DHS as epigenetic memory modules that support an accessible chromatin environment and thereby prime inducible genes for efficient re-activation when memory cells re-encounter Ag stimuli. In the absence of the priming elements, the lL-3/GM-CSF locus appears to remain inaccessible to transcription factors that are otherwise very efficient at recruiting the chromatin remodelers that create DHSs within the inducible enhancers.

Author details

¹Leeds Institute of Molecular Medicine, University of Leeds, Leeds LS9 7TF, UK. ²College of Medicine and Dentistry, Institute for Biomedical Research, University of Birmingham, Birmingham B152TT, UK.

Published: 18 March 2013

References

- Mirabella , et al: The human IL-3/GM-CSF locus is epigenetically silent in immature thymocytes and is progressively activated during T cell development. J Immunol 2010, 184:3043-3054.
- Baxter, et al: The inducible tissue-specific expression of the human IL-3/ GM-CSF locus is controlled by a complex array of developmentally regulated enhancers. J Immunol 2012, 189:4459-4469.

doi:10.1186/1756-8935-6-S1-O38

Cite this article as: James *et al.*: Prior epigenetic priming of cytokine genes in naive T cells is required for their subsequent activation by inducible enhancers. *Epigenetics & Chromatin* 2013 **6**(Suppl 1):O38.

Full list of author information is available at the end of the article



 $[\]overline{\ }^{1}\text{Leeds}$ Institute of Molecular Medicine, University of Leeds, Leeds LS9 7TF, UK